

# Recipe: Preparation of high performance chemical competent Cell

Pan M. CHU

September 14, 2020@ SIAT

Date Performed: August 7, 2019

## 1 Introduction

## 2 Protocol

- Streak glycerol stocks in plate with appropriate antibiotics;
- Pick one to three colonies into LB medium with appropriate antibiotics;
- Overnight culture dilute 100 fold in 50 mL SOB medium with appropriate antibiotics;
- Vigorously shaking before  $OD_{600}$  reach to about 0.6 (It is recommended having pre-culture in SOB medium for  $OD_{600}$  reach to 0.3.);
- Transform 50 mL medium into 50 mL centrifuge tube and stall it in ice for 10 min;
- 4 °C, 2500 ×g centrifuge 10 min;
- Discard supernatant, re-suspend pellet by 15 mL TB buffer (Table ??.)and stall it in ice for 10 min;
- 4 °C, 2500 ×g centrifuge 10 min;
- Discard supernatant and use 4 mL TB buffer re-suspend pellet, add 300  $\mu$ L DMSO finally(final concentration: 7% (v/v)).
- Add 5  $\mu$ L plasmid solution into 100  $\mu$ L fresh chemical competent cell.
- Stall on ice for 30 min.
- 42 °C heat-shocked for 45s and chilled on ice for 2 min.
- Add 900  $\mu$ L SOC medium and shake the culture vigorously.
- Coat 100  $\mu$ L medium on plate with appropriate antibiotics.

Component	Volume
ddH <sub>2</sub> O	12.5 mL
1 M KCl	4 mL
0.45M MnCl <sub>2</sub>	2.4 mL
0.5 M CaCl <sub>2</sub>	0.6 mL
0.5 M K-MES <sup>(1)</sup>	0.5 mL

Table 1: 20 mL K-MES Buffer Recipe

### Note:

- (1) Use KOH adjusting the K-MES solution to pH 6.3, store at −20°C for long term storage and split into aliquots avoiding repeated freezing and thawing.
- (2) For a 5 mL system. It can be performed *via* reducing volume of each reagent proportionally.